

WC 4



EASILY CONVERTS FROM A TRANSPORT CHAIR  
TO A SELF PROPELLED WHEELCHAIR.

OPTIONAL  
LEG RESTS  
AVAILABLE

SELF PROPELLED  
WHEELCHAIR



TRANSPORT  
CHAIR



aluminum transport chairs

To order, call toll free: 800-252-3054

drive  
MEDICAL DESIGN & MANUFACTURING

### Polly-Fly High Strength, Light Weight Wheelchair/Flyweight Transport Chair Combo

TR16 16" Seat Width, Silver Frame and Black Upholstery, 1/cs

TR18 18" Seat Width, Silver Frame and Black Upholstery, 1/cs

TR20 20" Seat Width, Silver Frame and Black Upholstery, 1/cs

LELR Swing-away Elevated Legrest, (Optional Front Riggings), 1/cs

#### FEATURES:

- Can be used as standard self propelled wheelchair or a transport chair all in one.
- Quick release 24" wheels can be removed to transition self propelled chair to a transport chair.
- All aluminum frame.
- Comes with two sets of wheel locks: One for use with wheelchair, one for use with transport chair.
- Fold down back.
- Nylon upholstery.
- Comes standard with swing away foot rests with composite footplates and heel loops.
- 8" casters
- Comes with anti-tippers.
- Limited Lifetime Warranty.

#### QUICK RELEASE ON 24" WHEEL.



#### SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/ Riggings	Product Weight	Carton Shipping Weight	Weight Limit
TR16	23.5"	12"	37"	19.5"	16"	16"	16"	42"	38 lbs.	43 lbs.	250 lbs.
TR18	24.5"	12"	37"	19.5"	16"	18"	16"	42"	38 lbs.	43 lbs.	250 lbs.
TR20	25.5"	13"	39.5"	19.5"	16"	20"	16"	41"	49 lbs.	56 lbs.	250 lbs.

Parkwood Oxygen & Medical Equipment  
7920 Congress Street Port Richey, FL 34668  
727-849-2577 or 800-252-3054